

Recommended Display Controllers

Status June 2016

Explanations

—	not recommended
•	recommended
•¹	additional DP-cable required



		Recommended Medical Display Controller						
		AMD			NVIDIA			
		FirePro 14.301.1010			Quadro 362.13			
		Driver Version ▶						
Model ▼		MD-W4100E	MD-W5100E	MD-W7100E	K-620E	K-2200E	M-2000E	MD-M4000E
Diagnostic Grayscale	MD211G3	•	•	•	•	•	•	•
	MD212G3	•	•	•	•	•	•	•
	MD211G5	—	•	•	—	•	•	•
Diagnostic Colour	MD21xC2	•	•	•	•	•	•	•
	MD21xC3	•	•	•	•	•	•	•
	MD242C2	•	•	•	•	•	•	•
	MD302C4	• ¹	• ¹	•	•	•	•	•
	MD302C6	—	•	•	—	•	•	•
	MD322C8 / -2	—	•	•	—	•	•	•
	MDC212C2	•	•	•	•	•	•	•
Review	MDview232	•	•	•	•	•	•	•
	MDview243	•	•	•	•	•	•	•
	MDview272	•	•	•	•	•	•	•
Form Factor		Half Length	Half Length	Full Length	Low Profile	Half Length	Half Length	Full Length
Display Connectors		4x mDP	4x DP	4x DP	1x DP; 1x DVI	2x DP; 1x DVI	4x DP	4x DP
Bus Connector		PCIe x16			PCIe x16			

Recommended Display Controllers

Status April 2015

Explanations	
—	not recommended
•	recommended
•¹	additional DP-cable required



		Recommended Medical Display Controller						
		AMD			Matrox		NVIDIA	
		FirePro 14.301.1010			Unified 2.10.1.6-a		Quadro 347.25	
Driver Version ►								
Model ▼		MD-4100E*	MD-W5100E *	MD-W7000E / W7100*	XENIA	XENIA PRO	MD-K600E / K620*	MD-K2000E / K2200*
Diagnostic Grayscale	MD211G3	●	●	●	●	●	●	●
	MD212G3	●	●	●	●	●	●	●
	MD211G5	—	●	●	—	●	—	●
Diagnostic Colour	MD21xC2	●	●	●	●	●	●	●
	MD21xC3	●	●	●	●	●	●	●
	MD242C2	●	●	●	●	●	●	●
	MD302C4	● ¹	● ¹	●	—	●	●	●
	MD302C6	—	●	●	—	●	—	●
	MD322C8	—	●	●	—	—	—	●
Review	MDview232	●	●	●	—	●	●	●
	MDview243	●	●	●	—	●	●	●
	MDview272	●	●	●	—	●	●	●
Form Factor		Half Length	Half Length	Full Length	Half Length	Half Length	Low Profile	Half Length
Display Connectors		4x mDP	4x DP	4x DP	3x DVI	3x DVI	1x DP; 1x DVI	2x DP; 1x DVI
Bus Connector		PCIe x16			PCIe x16		PCIe x16	

* Not supporting Windows XP